

Product Datasheet

Human C1QTNF1 ELISA Kit (Ready to Use) (orb550047)

Description	Human C1q And Tumor Necrosis Factor Related Protein 1 (C1QTNF1) Sandwich ELISA KitReady to Use: Pre- ...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9BXJ1
Sensitivity	0.57ng/mL
Entrez	114897
Expiration Date	Please enquire.