

Product Datasheet

Human CALPRO ELISA Kit (Ready to Use) (orb550086)

Description	Human Calprotectin (CALPRO) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conve...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P06702, P05109
Sensitivity	13.12pg/mL
Entrez	6279
Expiration Date	Please enquire.