

Product Datasheet

Human CA19-5 ELISA Kit (Ready to Use) (orb550100)

Description	Human Carbohydrate Antigen 19-5 (CA19-5) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.86ng/mL
Expiration Date	Please enquire.