

Product Datasheet

Human CEACAM6 ELISA Kit (Ready to Use) (orb550131)

Description	Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 6 (CEACAM6) Sandwich ELISA Kit Ready to...
Reactivity	Human
Range	23.4-1500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P40199
Sensitivity	9.68pg/mL
Entrez	4680
Expiration Date	Please enquire.