

## Product Datasheet

### Human CEACAM8 ELISA Kit (Ready to Use) (orb550133)

<b>Description</b>	Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 8 (CEACAM8) Sandwich ELISA Kit Ready to...
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<b>P31997</b>
<b>Sensitivity</b>	5.9pg/mL
<b>Entrez</b>	<b>1088</b>
<b>Expiration Date</b>	Please enquire.