

Product Datasheet

Human CASP14 ELISA Kit (Ready to Use) (orb550149)

Description	Human Caspase 14 (CASP14) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P31944
Sensitivity	0.106ng/mL
Entrez	23581
Expiration Date	Please enquire.