

Product Datasheet

Human Caveolin 1 ELISA Kit (Ready to Use) (orb550178)

Description	Human Caveolin 1 (Caveolin 1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	0.234-15ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q03135
Sensitivity	0.07ng/mL
Entrez	857
Expiration Date	Please enquire.