

Product Datasheet

Human CD27BP ELISA Kit (Ready to Use) (orb550184)

Description	Human CD27 Binding Protein (CD27BP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents f...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O15304
Sensitivity	0.068ng/mL
Entrez	10572
Expiration Date	Please enquire.