

Product Datasheet

Human CD5L ELISA Kit (Ready to Use) (orb550185)

Description	Human CD5 Antigen Like Protein (CD5L) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents...
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O43866
Sensitivity	28.6pg/mL
Entrez	922
Expiration Date	Please enquire.