

## Product Datasheet

### Human CDNF ELISA Kit (Ready to Use) (orb550197)

<b>Description</b>	Human Cerebral Dopamine Neurotrophic Factor (CDNF) Sandwich ELISA KitReady to Use: Pre-diluted detec...
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q49AH0</a>
<b>Sensitivity</b>	0.057ng/mL
<b>Entrez</b>	<a href="#">441549</a>
<b>Expiration Date</b>	Please enquire.