

## Product Datasheet

### APOL2 Immunizing peptide (orb55020)

<b>Description</b>	APOL2 Immunizing peptide
<b>Preservatives</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Storage</b>	Shipped at ambient temperature, store at - 20°C
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb55020 (APOL2 antibody)
<b>Protein Sequence</b>	TKIHMLQPGQD
<b>Expiration Date</b>	12 months from date of receipt.