

Product Datasheet

Human CXCR4 ELISA Kit (Ready to Use) (orb550215)

Description	Human Chemokine C-X-C-Motif Receptor 4 (CXCR4) Sandwich ELISA KitReady to Use: Pre-diluted detection...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P61073
Sensitivity	0.060ng/mL
Entrez	7852
Expiration Date	Please enquire.