

Product Datasheet

Human CCKAR ELISA Kit (Ready to Use) (orb550227)

Description	Human Cholecystokinin A Receptor (CCKAR) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P32238
Sensitivity	0.084ng/mL
Entrez	886
Expiration Date	Please enquire.