

Product Datasheet

Human CHRNa3 ELISA Kit (Ready to Use) (orb550234)

| | |
|--------------------------|---|
| Description | Human Cholinergic Receptor, Nicotinic, Alpha 3 (CHRNa3) Sandwich ELISA KitReady to Use: Pre-diluted ... |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids. |
| Sample Types | tissue homogenates and other biological fluids. |
| Uniprot ID | P32297 |
| Sensitivity | 0.107ng/mL |
| Entrez | 1136 |
| Expiration Date | Please enquire. |