

Product Datasheet

Human CHGA ELISA Kit (Ready to Use) (orb550246)

Description	Human Chromogranin A (CHGA) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conve...
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P10645
Sensitivity	23.5pg/mL
Entrez	1113
Expiration Date	Please enquire.