

Product Datasheet

Human CLDN3 ELISA Kit (Ready to Use) (orb550255)

| | |
|--------------------------|--------------------------------------------------------------------------------------------------------------------|
| Description | Human Claudin 3 (CLDN3) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien... |
| Reactivity | Human |
| Range | 0.312-20ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | O15551 |
| Sensitivity | 0.125ng/mL |
| Entrez | 1365 |
| Expiration Date | Please enquire. |