

Product Datasheet

Human CD83 ELISA Kit (Ready to Use) (orb550278)

Description	Human Cluster Of Differentiation 83 (CD83) Sandwich ELISA Kit Ready to Use: Pre-diluted detection rea...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, cell lysates or other biological fluids.
Sample Types	serum, plasma, cell lysates and other biological fluids.
Uniprot ID	Q01151
Sensitivity	13.2pg/mL
Entrez	9308
Expiration Date	Please enquire.