

Product Datasheet

Human COL10 ELISA Kit (Ready to Use) (orb550312)

Description	Human Collagen Type X (COL10) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q03692
Sensitivity	0.141ng/mL
Entrez	1300
Expiration Date	Please enquire.