

Product Datasheet

Human CFH ELISA Kit (Ready to Use) (orb550347)

Description	Human Complement Factor H (CFH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c...
Reactivity	Human
Range	2.34-150ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P08603
Sensitivity	0.85ng/mL
Entrez	3075
Expiration Date	Please enquire.