

Product Datasheet

Human CPLX2 ELISA Kit (Ready to Use) (orb550355)

Description	Human Complexin 2 (CPLX2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6PUV4
Sensitivity	0.084ng/mL
Entrez	10814
Expiration Date	Please enquire.