

Product Datasheet

Human CORT ELISA Kit (Ready to Use) (orb550375)

| | |
|--------------------------|---|
| Description | Human Cortistatin (CORT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie... |
| Reactivity | Human |
| Range | 1.56-100ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | 000230 |
| Sensitivity | 0.66ng/mL |
| Entrez | 1325 |
| Expiration Date | Please enquire. |