

Product Datasheet

Human CKB ELISA Kit (Ready to Use) (orb550379)

Description	Human Creatine Kinase, Brain (CKB) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents fo...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P12277
Sensitivity	0.58ng/mL
Entrez	1152
Expiration Date	Please enquire.