

Product Datasheet

Human CRT ELISA Kit (Ready to Use) (orb550382)

Description	Human Creatine Transporter (CRT) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for ...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P48029
Sensitivity	0.124ng/mL
Entrez	6535
Expiration Date	Please enquire.