

## Product Datasheet

# Human CTXIII ELISA Kit (Ready to Use) (orb550388)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human Cross Linked C-Telopeptide Of Type III Collagen (CTXIII) Sandwich ELISA KitReady to Use: Pre-d... |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 0.156-10ng/mL   |
| <b>Form/Appearance</b>   | Ready to Use  |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.          |
| <b>Sample Types</b>      | serum, plasma and other biological fluids.  |
| <b>Sensitivity</b>       | 0.043ng/mL  |
| <b>Expiration Date</b>   | Please enquire.   |