

Product Datasheet

Human CYPB ELISA Kit (Ready to Use) (orb550412)

Description	Human Cyclophilin B (CYPB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P23284
Sensitivity	32pg/mL
Entrez	5479
Expiration Date	Please enquire.