

Product Datasheet

Human CRISPLD2 ELISA Kit (Ready to Use) (orb550430)

Description	Human Cysteine-Rich Secretory Protein LCCL Domain-Containing 2 (CRISPLD2) Sandwich ELISA KitReady to...
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	1.12ng/mL
Expiration Date	Please enquire.