

Product Datasheet

Human CYP2B6 ELISA Kit (Ready to Use) (orb550449)

Description	Human Cytochrome P450 2B6 (CYP2B6) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents fo...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P20813
Sensitivity	0.31ng/mL
Entrez	1555
Expiration Date	Please enquire.