

Product Datasheet

Human DMBT1 ELISA Kit (Ready to Use) (orb550490)

Description	Human Deleted In Malignant Brain Tumors 1 (DMBT1) Sandwich ELISA Kit Ready to Use: Pre-diluted detect...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, lung lavage fluid or other biological fluids.
Sample Types	tissue homogenates, lung lavage fluid and other biological fluids.
Uniprot ID	Q9UGM3
Sensitivity	0.058ng/mL
Entrez	1755
Expiration Date	Please enquire.