

Product Datasheet

Human DYSF ELISA Kit (Ready to Use) (orb550581)

Description	Human Dysferlin (DYSF) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	075923
Sensitivity	0.141ng/mL
Entrez	8291
Expiration Date	Please enquire.