

Product Datasheet

Human EDIL3 ELISA Kit (Ready to Use) (orb550600)

Description	Human EGF Like Repeats And Discoidin I Like Domains Protein 3 (EDIL3) Sandwich ELISA KitReady to Use...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O43854
Sensitivity	0.131ng/mL
Entrez	10085
Expiration Date	Please enquire.