

Product Datasheet

Human EDN2 ELISA Kit (Ready to Use) (orb550628)

Description	Human Endothelin 2 (EDN2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P20800
Sensitivity	11.6pg/mL
Entrez	1907
Expiration Date	Please enquire.