

Product Datasheet

Human ECE1 ELISA Kit (Ready to Use) (orb550629)

Description	Human Endothelin Converting Enzyme 1 (ECE1) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P42892
Sensitivity	11.3pg/mL
Entrez	1889
Expiration Date	Please enquire.