

Product Datasheet

Human EPX ELISA Kit (Ready to Use) (orb550641)

| | |
|--------------------------|---|
| Description | Human Eosinophil Peroxidase (EPX) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for... |
| Reactivity | Human |
| Range | 1.25-80ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | P11678 |
| Sensitivity | 0.51ng/mL |
| Entrez | 8288 |
| Expiration Date | Please enquire. |