

Product Datasheet

Human EPHA10 ELISA Kit (Ready to Use) (orb550648)

Description	Human Ephrin Type A Receptor 10 (EPHA10) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q5JZY3
Sensitivity	0.109ng/mL
Entrez	284656
Expiration Date	Please enquire.