

Product Datasheet

Human EGFR2 ELISA Kit (Ready to Use) (orb550651)

Description	Human Epidermal Growth Factor Receptor 2 (EGFR2) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P04626
Sensitivity	0.25ng/mL
Entrez	2064
Expiration Date	Please enquire.