

Product Datasheet

Human FcgBP ELISA Kit (Ready to Use) (orb550717)

Description	Human Fc Fragment Of IgG Binding Protein (FcgBP) Sandwich ELISA KitReady to Use: Pre-diluted detecti...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9Y6R7
Sensitivity	0.116ng/mL
Entrez	8857
Expiration Date	Please enquire.