

Product Datasheet

Human FGF13 ELISA Kit (Ready to Use) (orb550738)

Description	Human Fibroblast Growth Factor 13 (FGF13) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q92913
Sensitivity	5.4pg/mL
Entrez	2258
Expiration Date	Please enquire.