

## Product Datasheet

# Human FGFR2 ELISA Kit (Ready to Use) (orb550750)

<b>Description</b>	Human Fibroblast Growth Factor Receptor 2 (FGFR2) Sandwich ELISA KitReady to Use: Pre-diluted detect...
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<b>P21802</b>
<b>Sensitivity</b>	12.2pg/mL
<b>Entrez</b>	<b>2263</b>
<b>Expiration Date</b>	Please enquire.