

Product Datasheet

Human FBLN4 ELISA Kit (Ready to Use) (orb550761)

Description	Human Fibulin 4 (FBLN4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Human
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	095967
Sensitivity	1.48ng/mL
Entrez	30008
Expiration Date	Please enquire.