

Product Datasheet

Human FLII ELISA Kit (Ready to Use) (orb550779)

Description	Human Flightless I Homolog (FLII) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q13045
Sensitivity	29.2pg/mL
Entrez	2314
Expiration Date	Please enquire.