

Product Datasheet

Human FLOT2 ELISA Kit (Ready to Use) (orb550780)

Description	Human Flotillin 2 (FLOT2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q14254
Sensitivity	0.216ng/mL
Entrez	2319
Expiration Date	Please enquire.