

Product Datasheet

Human FOLR2 ELISA Kit (Ready to Use) (orb550785)

| | |
|--------------------------|--|
| Description | Human Folate Receptor 2, Adult (FOLR2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent... |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | P14207 |
| Sensitivity | 0.064ng/mL |
| Entrez | 2350 |
| Expiration Date | Please enquire. |