

Product Datasheet

Human gGT5 ELISA Kit (Ready to Use) (orb550837)

Description	Human Gamma-Glutamyltransferase 5 (gGT5) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P36269
Sensitivity	0.55ng/mL
Entrez	2687
Expiration Date	Please enquire.