

Product Datasheet

Human GKN2 ELISA Kit (Ready to Use) (orb550844)

Description	Human Gastrokine 2 (GKN2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q86XP6
Sensitivity	6.1pg/mL
Entrez	200504
Expiration Date	Please enquire.