

Product Datasheet

Human GFRa1 ELISA Kit (Ready to Use) (orb550854)

Description	Human Glial Cell Line Derived Neurotrophic Factor Receptor Alpha 1 (GFRa1) Sandwich ELISA KitReady t...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P56159
Sensitivity	0.119ng/mL
Entrez	2674
Expiration Date	Please enquire.