

Product Datasheet

Human GZMM ELISA Kit (Ready to Use) (orb550948)

Description	Human Granzyme M (GZMM) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P51124
Sensitivity	5.9pg/mL
Entrez	3004
Expiration Date	Please enquire.