

Product Datasheet

Human HSPB1 ELISA Kit (Ready to Use) (orb550979)

Description	Human Heat Shock 27kDa Protein 1 (HSPB1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P04792
Sensitivity	0.32ng/mL
Entrez	3315
Expiration Date	Please enquire.