

Product Datasheet

Human HSPBP1 ELISA Kit (Ready to Use) (orb550982)

Description	Human Heat Shock 70kDa Binding Protein 1 (HSPBP1) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9NZL4
Sensitivity	0.105ng/mL
Entrez	23640
Expiration Date	Please enquire.