

Product Datasheet

Human HSPA8 ELISA Kit (Ready to Use) (orb550988)

Description	Human Heat Shock 70kDa Protein 8 (HSPA8) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P11142
Sensitivity	0.134ng/mL
Entrez	3312
Expiration Date	Please enquire.