

Product Datasheet

Human HSPb2 ELISA Kit (Ready to Use) (orb550994)

Description	Human Heat Shock Protein Beta 2 (HSPb2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagen...
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q16082
Sensitivity	0.244ng/mL
Entrez	3316
Expiration Date	Please enquire.